

FRONT VIEW

MEDICAL GAS OUTLET PANEL:

- PANEL SIZE VARIES DEPENDING ON NUMBER OF OUTLETS
- CONFIGURATION OF OUTLETS NEEDS TO BE ADVISED

WHEN OUTLET PANEL IN CUPBOARD, THE FOLLOWING IS REQUIRED:

- MINIMUM 200MM CLEARANCE ABOVE
- MINIMUM 300MM CLEARANCE BELOW
- NEED TO ALLOW FOR CATCH POT TO BE OUTSIDE THE CUPBOARD

FOR APPROVAL



NOTE ADDED 29-02-16 17-02-16 FOR APPROVAL Description Checked Project SGI STANDARD FITTINGS SGI MGOP-3G.DWG Drawn DGL 17/02/2016 MEDICAL GAS OUTLET PANEL Ch'ked Scale Drawing Number 3 GANG MGOP-3G В

Property of Special Gases Enterprises Pty Ltd